



**Comments for pZeoSV2/lacZ**  
**6661 nucleotides**

- SV40 promoter: bases 10-356
- T7 promoter/priming site: 401-420
- LacZ ORF: bases 458-3626
- BGH polyA: bases 3706-3919
- pUC origin: bases 3959-4683
- f1 origin: bases 4926-5436
- CMV promoter: bases 5450-6076
- EM-7 promoter: bases 6077-6143
- Zeocin resistance gene ORF: bases 6144-6518
- SV40 polyA: bases 6526-6661