


Thermo Scientific Dionex Integrion HPIC System

Remote Services

Product Spotlight

The Thermo Scientific™ Dionex™ Integrion™ HPIC System brings the power of Unity™ Remote Services to your ion chromatography system, increasing your productivity by maximizing instrument uptime and reducing costs.

Unity Remote Services software enables remote support, monitoring of optimal instrument operation, and remote detection and diagnosis of issues. Hardware performance metrics are monitored over the internet to expand instrument support beyond just care and maintenance.



Maximizing instrument uptime

- Allows Unity Lab Services to proactively identify and respond to existing and potential problems quickly and easily
- Avoids the need for your personnel to spend time with trial-and-error troubleshooting, reducing costs

Reducing costs

- Proactively assesses upcoming maintenance needs and optimizes timing to fit with your lab's flow
- Reduces the need for unscheduled service visits
- Avoids the loss of your samples and critical data

Complete data security

- Complies with multiple regulatory requirements (HIPAA, CLSI Auto 9A, ISO 27001)
- Uses the latest HTTPS, PKI, and up to 256-bit AES encryption technology to securely transfer data from your lab
- Firewall friendly; uses standard HTTPS on port 443
- Allows for access to instrument system diagnostic information only

For more information, visit thermoscientific.com

©2015 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Unity Lab Services

A Thermo Fisher Scientific Brand

Thermo
SCIENTIFIC

A Thermo Fisher Scientific Brand