



Precision is a perfect fit

Lower the cost of consumables and reduce waste with one plate that fits two block formats: 96- and 32-well



The new Applied Biosystems™ MicroAmp™ TriFlex 3 x 32-Well PCR Reaction Plate is designed for optimal performance with the Applied Biosystems™ ProFlex™ 3 x 32-Well PCR System and all Applied Biosystems™ standard 96-well blocks. Just tear off what you need and get on with your work.

- One plate that can be used as separated, individual 32-well plates or as a full 96-well plate
- Easy-to-tear dual side tabs for simple plate separation
- Alphanumeric labeling on each plate segment for easy well identification
- Applied Biosystems™ MicroAmp™ 32-Well Clear Adhesive Film for evaporation-free plate sealing

Designed specifically for Applied Biosystems™ thermal cyclers to enable optimal PCR results.

appliedbiosystems
by Thermo Fisher Scientific

Engineer Approved

All Applied Biosystems thermal cyclers and PCR plastics are validated to help ensure optimal PCR performance, and are “Engineer Approved”—designed and built for reliability and precision to enable optimal PCR results.

Get the most from your low- and standard-throughput PCR instruments by choosing the MicroAmp TriFlex 3 x 32-Well PCR Reaction Plate for 96 or 32 wells—one plate fits all



Applied Biosystems™ SimpliAmp™
Thermal Cycler



Applied Biosystems™ ProFlex™
3 x 32-well PCR System

Designed for use with
the ProFlex 3 x 32-Well
PCR System



Applied Biosystems™ ProFlex™
96-well PCR System

Compatible with all Applied
Biosystems standard 96-well
block thermal cyclers



MicroAmp 32-Well Clear
Adhesive Film for evaporation-
free plate sealing

Instrument compatibility

Product	Cat. No.	3 x 32-well		96-well	
		ProFlex	Veriti, ProFlex, SimpliAmp	2720	9700
96-well 0.2 mL plates					
MicroAmp TriFlex 3 x 32-Well PCR Reaction Plate	A32810, A32811	•	•	•	•
Seals and covers					
MicroAmp 32-Well Clear Adhesive Film	A32812	•	•	•	•

Ordering information

Product	Quantity	Cat. No.
MicroAmp TriFlex 3 x 32-Well PCR Reaction Plate	5 plates/pack	A32810
MicroAmp TriFlex 3 x 32-Well PCR Reaction Plate	20 plates/pack	A32811
MicroAmp 32-Well Clear Adhesive Film	100 films/pack	A32812

Reduce cost and waste at thermofisher.com/3x32plate

ThermoFisher
SCIENTIFIC